

Large Capacity Rotors

# Thermo Scientific 6 x 2000 mL Swinging Bucket Rotor

6 x 2000 mL



## ROTOR SPECIFICATIONS

<b>Type</b>	Swinging Bucket
<b>Material</b>	Stainless Steel Rotor Body with Nickel-Coating; Aluminum Windshield and Buckets <sup>1</sup> ; Polypropylene Bottles
<b>Capacity</b>	6 x 2000 mL
<b>Maximum Speed</b>	4700 rpm
<b>Bucket Cavity Dimensions</b>	100 x 185 x 150 mm
<b>Angle</b>	90°
<b>K Factor</b>	7909
<b>Net Weight</b>	
Bucket	1.8 kg / 3.9 lbs
Windshield	5.4 kg / 11.9 lbs
Lid	1.7 kg / 3.7 lbs
Rotor Body	18.5 kg / 40.7 lbs
<b>Warranty</b>	7 years

	RCF	RADIUS
<b>Maximum</b>	7187 x g	291 mm
<b>Average</b>	5396 x g	219 mm
<b>Minimum</b>	3605 x g	146 mm

## ROTOR AND ACCESSORIES ORDERING INFORMATION

		THERMO SCIENTIFIC SORVALL BIOS 16
<b>6 x 2000 mL Rotor</b>	Max. Speed	4700 rpm
	Max. RCF	7187 x g

### 6 x 2000 mL Rotor Packages include:

Description	Cat. No.
6 x 2000 mL Swinging Bucket Rotor, Oval Buckets (qty. 6) and 2000 mL Bottles (qty. 6)	75003962
Rotor Kit	70900140
Rotor Locking Pin	75003856
Rotor Locking Tool	11246
Bolt Grease	75003789
Corrosion Protection Oil	70009824
Cleaning Pad (qty. 6)	66309
Information Card	50149709
6 x 2000 mL Swinging Bucket Rotor, Oval Buckets with Dura-Coat (qty. 6) and 2000 mL Bottles (qty. 6)	75004162
Rotor Kit	70900140
Rotor Locking Pin	75003856
Rotor Locking Tool	11246
Bolt Grease	75003789
Corrosion Protection Oil	70009824
Cleaning Pad (qty. 6)	66309
Information Card	50149709

supplyLAB

www.supplylab.pt  
geral@supplylab.pt

Cacém Park - Edifício 9  
Estrada de Paço de Arcos nº88  
2739-512 Aqualva Cacém  
T +(351) 21 4278700  
F +(351) 21 4278709

**Thermo**  
SCIENTIFIC

**ROTOR AND ACCESSORIES ORDERING INFORMATION**
**6 x 2000 mL Rotor Accessories**

Description	Cat. No.
Oval Bucket (qty. 2)	75003964

**TUBE AND BOTTLE ORDERING INFORMATION**

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
2000.0	–	Oval Bottle	75003870	2	6	4,700	7187	–	Closure	Included	2	PP Sealing Lid

**supplyLAB**

www.supplylab.pt  
geral@supplylab.pt

Cacém Park - Edifício 9  
Estrada de Paço de Arcos nº88  
2739-512 Agualva Cacém  
T +(351) 21 4278700  
F +(351) 21 4278709

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**For more information, please visit [www.thermofisher.com/rotors](http://www.thermofisher.com/rotors).**