



TFT-80.2¹

Titanium Rotor
12-Place, 2 ml
Tube Dimensions: 8 x 51 mm
Fixed-Angle: 25°
Maximum Speed: 80,000 rpm²
K Factor: 22.6
Rotor Mass: 2.5 kg

	RCF (g)	R (cm)
Max.	429,643	6.01
Avg.	335,994	4.700
Min.	242,344	3.39

Basic Catalog No. 54456	Complete Catalog No. 54457	Accessories
Includes: TFT-80.2 Rotor Rotor Cover 54463 Locking Nut 54464 Rotor Stand 51942 Vacuum Grease 65937 Lubricant 61556 Rotor Handle 11559 Overspeed Decal, 80,000 rpm, Extra 11628 Rotor Cover O-ring, Viton®, Extra (5 supplied) 54465 Locking Nut O-ring, Viton®, Extra (5 supplied) 54466 ULTRASPEED Rotor Log Book 52384 Instruction Manual 54462	Includes Basic Catalog No. 54456 Plus: Tubes, Polyallomer, 2 ml (4 pkgs, 50/pkg) 54460 Sealing Cap Assembly, Multipiece, Selfseal Caps (1per package) 54459	Drive Adapter Kit for use in some RC-ULTRA Centrifuges 54381 Tube Removal Tool 11538

TUBES

CAPACITY (ml)		DESCRIPTION	MAX SPEED ² (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	CAT. NO. (Qty/Pkg)	NO. OF PLACES	
2.0	1.3	Polyallomer Thin-Walled Tube	80,000	8 x 51	54460 (50)	54459 (1)	Multipiece, Aluminum Selfseal	—	—	Requires Tube Removal Tool 11538 .

⁽¹⁾ May require Drive Adapter (54381) for RC-ULTRA centrifuges.

⁽²⁾ Non-precipitating solutions up to 1.2 g/ml can be run in this rotor without a reduction in rotor speed.

All trademarks are the property of their respective owners.

supplyLAB

www.supplylab.pt
geral@supplylab.pt

Cacém Park - Edifício 9
Estrada de Paço de Arcos nº88
2739-512 Agualva Cacém
T +(351) 21 4278700
F +(351) 21 4278709