



TSB Series blood bank refrigerators

Enhanced functionality for secure
blood storage

Security and quality for environmentally sustainable storage of blood and blood components

Thermo Scientific™ TSB Series blood bank refrigerators are designed to meet the strict requirements for storage of blood and blood-related products in full compliance with 93/42/CE—class IIa regulation.



Security

- Digital touchscreen control
- Onboard datalogging tracking temperatures and alarms
- Electronic lock with password protection
- Power failure alarm with battery backup
- Ballasted probe to simulate product temperature




Quality

- 12-month warranty
- 93/42/CE—class IIa regulation
- White painted interior and precoated steel exterior
- Robust 75 mm thick insulation



Sustainability

- F-gas compliant hydrocarbon green refrigerants
- Low energy consumption with real-time energy monitoring
- Low heat output



TSB Series blood bank refrigerators offer:

- A wide range of sizes from 140 L undercounter to 1,500 L double door
- Touchscreen digital controller with onboard datalogging
- Whisper-quiet operation of <50 dBA
- Triple-pane glass doors
- White painted interior and exterior
- Temperature, door ajar, and power failure (with battery backup) alarms
- Low energy consumption with real-time energy monitoring
- Low heat output
- Rear access port
- Electronic lock with password protection

Ordering information

Capacity (L)	Doors	Drawers	Dimensions (W x D x H)	Voltage	Plug type	Packaging	Cat. No.
140	1 glass	3	520 x 635 x 1,000 mm	230 V/50 Hz	EU	Box	TSB140REGBV
						Crate	TSB140REGCV
					UK	Box	TSB140REGBW
						Crate	TSB140REGCW
250		5	600 x 655 x 1,635 mm		EU	Box	TSB250REGBV
						Crate	TSB250REGCV
					UK	Box	TSB250REGBW
						Crate	TSB250REGCW
400		6	600 x 655 x 1,955 mm		EU	Box	TSB400REGBV
						Crate	TSB400REGCV
					UK	Box	TSB400REGBW
						Crate	TSB400REGCW
700	7	750 x 840 x 2,100 mm	EU	Box	TSB700REGBV		
				Crate	TSB700REGCV		
			UK	Box	TSB700REGBW		
				Crate	TSB700REGCW		
1,500	2 glass	14	1440 x 825 x 2,070 mm	EU	Box	TSB1500REGBV	
					Crate	TSB1500REGCV	
				UK	Box	TSB1500REGBW	
					Crate	TSB1500REGCW	

Find out more at thermofisher.com/cold